




## Start File




1. Click "Next"




2. If you agree with the Software License Agreement please select "I Agree", then click "Next"




3. If you do not know what are these fields for please do not change anything just click "Next"



4. Instead of dentalsoftware, please enter your password which you have created during the installation of MySQLServer. Everything else should look like in the image. Click "Next" .




5. You have started the installation



6. When the installation finishes click "Next"


Icon should appear on desktop or in start menu , please click on the icon so you can start the program.

When you first start the software,it should look like the image below.




Depending on whether you want to buy or test our software, click “Generate Trial License Request” or “Generate Full License Request”

Fill the field “E-Mail” and “Phone” .




If you would like to send us your request via e-mail, click “Yes” , or if you would like to give us a call click “No”



When you generate a request a serial number should appear like on the image above please write that number down and either send it by mail or give us a call. Close the program!

After we send you the licence start the program again.



Enter Licence number and then click (in this case) Register Full Licence



Please click on “Test connection”




If you get a message “Test connection succeeded” , please click on “Create database”




If you get a message “Successfully created database!” , click “Close”.

Start your application again



Default username: user, password: pass

Enter that in the fields and press “Ok”



If an image appears as above you have installed the software successfully.