

Dental Software – Dental Office Assistance.


Thank you for downloading our software.

Please read the instructions very carefully ,if you do not follow the instructions carefully, the program will not work.


First start MYSQL Server.




1. Click “Next”




2. Select option “Complete” and then click “Next”.




3. Click "Install"




4. Window as above should appear




5. Click "Next"




6. Check if “Configure the MySQL Server now” has been checked, then click Finish




7. Click “Next”




8. Click "Install"




9. Check if the "Developer Machine" is checked then click "Next"




10. Check if “Multifunctional Database” has been checked, then click “Next”




11. Your window should look like the image above, please do not make any changes, click “Next”.




12. Check if "Decision Support(DSS)/OLAP" has been checked ,then click "Next"




13. All fields need to be checked !




14. Check if “Standard Character Set” has been checked, then click “Next”




15. Make sure that you check all the fields , if you do not check everything the program will not work.




16. Enter the password and then retype it again. Please remember your password.



17. When you enter your password, click "Next".



18. Then click Execute.



19. When execution is complete and all the fields have been checked automatically click "Finish".

My Sql Server has been installed